



FRINGS® & INOX Instruments

Instrument Management Guide

- *All stainless-steel instruments must be cleaned and thoroughly dried prior to sterilization.*
- *We recommend using only a non-corrosive detergent.*
- *Abrasive brushes or materials should not be used to clean instruments.*

STEP ONE:

Following clinical use, all instruments must be thoroughly rinsed with water to remove all residue prior to instrument processing.

STEP TWO:

Instruments should be placed into an ultrasonic bath for a period of 20-30 minutes according to manufacturer guidelines to ensure complete instrument cleaning. NOTE: Ultrasonic Solution should be changed daily or when the solution appears to be visibly soiled.

STEP THREE:

All instruments must be rinsed with water and completely dried before autoclaving. All stainless-steel instruments should be placed in separate sterilization pouches and not to be mixed with any carbide instruments during instrument processing.

NOTE: Only demineralized water should be used in autoclave equipment according to autoclave manufacturer instructions.

FRINGS® Maintenance and Care

All **FRINGS®** must be lubricated twice a month to ensure optimal performance. We recommend a surgical milk bath every two weeks as part of your instrument management process. Following the surgical milk bath, lubricate one drop of surgical instrument oil inside the hinge.